

## Isoalantolactone

Catalog No: #S2829


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Isoalantolactone
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	232.31
Formulation	C15H20O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:C[C@]12CCCC(=C)[C@@H]1C[C@H]3[C@@H](C2)OC(=O)C3=CPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Transparent Needle CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Di W, et al. Oncol Rep. 2014 Oct;32(4):1585-93.

Note: This product is for in vitro research use only