

## Glycocholic acid

Catalog No: #S2831


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Glycocholic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	465.64
Formulation	C26H43NO6
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OtherSMILES:C[C@H](CCC(=O)NCC(=O)O)[C@H]1CC[C@@H]2[C@@]1([C@H](C[C@H]3[C@H]2[C@@H](C[C@H]4[C@@]3(CC[C@H](C4)O)C)O)O)CPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: SolubleAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Zhang A, et al. Clin Chim Acta. 2013 Mar 15;418:86-90.

Note: This product is for in vitro research use only