

S-Isocorydine(+)

Catalog No: #S2832


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	S-Isocorydine(+)
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	341.41
Formulation	C ₂₀ H ₂₃ NO ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: CN1CCC2=CC(=C(C3=C2[C@@H]1CC4=C3C(=C(C=C4)OC)O)OC)OC
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Sun H, et al. PLoS One. 2012;7(5):e36808.

Note: This product is for in vitro research use only