

Dehydrocostus Lactone

Catalog No: #S2833


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dehydrocostus Lactone
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	230.3
Formulation	C15H18O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:C=C1CC[C@@H]2[C@@H]([C@@H]3[C@H]1CCC3=C)OC(=O)C2=CPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Dong GZ, et al. PhytOthers Res. 2015 May;29(5):680-6.

Note: This product is for in vitro research use only