

## Nobiletin

Catalog No: #S2834


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Nobiletin
Brief Description	Inhibitors
Purification	99.00%
Target Name	MMP inhibitor
Calculated MW	402.4
Formulation	C <sub>21</sub> H <sub>22</sub> O <sub>8</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Other SMILES: COC1=C(OC)C=C(C=C1)C1=CC(=O)C2=C(O1)C(OC)=C(OC)C(OC)=C2OC Pathways: Proteases/Proteasome Receptor: MMP Bolin  
 g pt: Melting pt: Solubility: Ethanol: 3 mg/mL (7.45 mM); DMSO: 81 mg/mL (201.29 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Murakami A, et al. Y Res, 2000, 60(18), 5059-5066.

Note: This product is for in vitro research use only