

Morin

Catalog No: #S2835


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Morin
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others inhibitor
Calculated MW	302.23
Formulation	C15H10O7
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:OC1=CC(O)=C(C(C(O)=C(C2=C(O)C=C(O)C=C2)O3)=O)C3=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Tian, WX. Current Medicinal Chemistry, 2006, 13 (8): 967C977.

Note: This product is for in vitro research use only