

Isorhamnetin

Catalog No: #S2836


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Isorhamnetin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	316.27
Formulation	C16H12O7
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: COC1=C(O)C=CC(=C1)C1=C(O)C(=O)C2=C(O)C=C(O)C=C2O1
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt:
 Solubility: DMSO: 10 mM
 Appearance: Yellow Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Han X, et al. *Biomol Ther (Seoul)*. 2015 Jul;23(4):357-66.

Note: This product is for in vitro research use only