

Naringenin


Catalog No: #S2838

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Naringenin
Brief Description	Inhibitors
Purification	99.00%
Target Name	P450 agonist
Calculated MW	272.26
Formulation	C15H12O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: C1[C@H](OC2=CC(=CC(=C2C1=O)O)O)C3=CC=C(C=C3)O
 Pathways: Metabolism
 Receptor: CYP1A2
 Boiling pt: Melting pt: Solubility: Ethanol: 5 mg/mL (18.36 mM); DMSO: 54 mg/mL (198.34 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Fuhr U, et al. Br J Clin Pharmacol, 1993, 35(4):431-436.

Note: This product is for in vitro research use only