

Orsellinic acid

Catalog No: #S2839


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Orsellinic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	168.15
Formulation	C ₈ H ₈ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CC1=CC(O)=CC(O)=C1C(=O)O
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: Gray Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Kang Ying, et al. Organic Letters, 2003, 5 (23): 4481C4484.

Note: This product is for in vitro research use only