

Phloracetophenone


Catalog No: #S2840

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phloracetophenone
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	168.15
Formulation	C ₈ H ₈ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CC(=O)C1=C(O)C=C(O)C=C1O
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Lai YH, et al. Rapid Commun Mass Spectrom. 2014 Aug 15;28(15):1716-22.

Note: This product is for in vitro research use only