

Tangeretin

Catalog No: #S2842


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tangeretin
Brief Description	Inhibitors
Purification	99.70%
Target Name	Others inhibitor
Calculated MW	372.38
Formulation	C ₂₀ H ₂₀ O ₇
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: COC1=CC=C(C=C1)C1=CC(=O)C2=C(O1)C(OC)=C(OC)C(OC)=C2OC
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: DMSO: 8 mg/mL (21.48 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

- Hirano T, et al. Br J Y, 1995, 72(6), 1380-1388.

Note: This product is for in vitro research use only