

## Hyperoside

Catalog No: #S2844


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Hyperoside
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	464.38
Formulation	C <sub>21</sub> H <sub>20</sub> O <sub>12</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Other SMILES: C1=CC(=C(C=C1)C2=C(C(=O)C3=C(C=C(C=C3O2)O)O)O)[C@H]4[C@@H]([C@H]([C@H]([C@H](O4)CO)O)O)O)O Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Zhang N, et al. PLoS One. 2014 Jul 1;9(7):e98973.

Note: This product is for in vitro research use only