

Indigo


Catalog No: #S2847

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Indigo
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	262.26
Formulation	C16H10O2N2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: C1=CC=C2C(=C1)C(=O)/C(=C\3/C(=O)C4=CC=CC=C4N3)/N2
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: Soluble
Appearance: Dark Blue Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Margetis K, et al. J Neurosurg. 2015 Jan;122(1):40-8.

Note: This product is for in vitro research use only