

Osthole

Catalog No: #S2848


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Osthole
Brief Description	Inhibitors
Purification	99.00%
Target Name	Histamine Receptor agonist
Calculated MW	244.29
Formulation	C15H16O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:COC1=C(CC=C(C)C)C2=C(C=CC(=O)O2)C=C1Pathways:GPCR/G Protein; Immunology/Inflammation;
NeuroscienceReceptor:H1 receptorBoiling pt:Melting pt:Solubility:Ethanol: 49 mg/mL (200.58 mM); DMSO: 49 mg/mL (200.58 mM)Appearance:White
PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be
stored below -20 °C for several month.

References

1. Ko FN, et al. Thromb Haemost, 1989, 62(3), 996-999.

Note: This product is for in vitro research use only