

Proanthocyanidins

Catalog No: #S2849


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Proanthocyanidins
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others
Calculated MW	594.52
Formulation	C30H26O13
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:OC1C(O)C(OC2CC3=C(O)C=C(O)C=C3OC2C2=CC(O)=C(O)C=C2)(OC2=CC(O)=CC(O)=C12)C1=CC(O)=C(O)C=C1Pathways:
OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Red PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Ma J, et al. Nan Fang Yi Ke Da Xue Xue Bao. 2015 Aug;35(8):1110-5. Chinese.

Note: This product is for in vitro research use only