

Citropten

Catalog No: #S2852


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Citropten
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	206.2
Formulation	C11H10O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:COC1=CC(OC)=C2C=CC(=O)OC2=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Zdybel M, et al. Acta Pol Pharm. 2014 Nov-Dec;71(6):1066-72.

Note: This product is for in vitro research use only