

Forsythin

Catalog No: #S2854


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Forsythin
Brief Description	Inhibitors
Purification	99%+
Target Name	p38 MAPK inhibitor
Calculated MW	534.56
Formulation	C27H34O11
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:COC1=C(OC)C=C(C=C1)C1OCC2C1COC2C1=C(OC)C=C(OC2OC(CO)C(O)C(O)C2O)C=C1Pathways:MAPKReceptor:p38

MAPKBoling pt:Melting pt:Solubility:DMSO: SolubleAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Pan X, et al. Inflamm Res. 2014 Jul;63(7):597-608.

Note: This product is for in vitro research use only