

Epicatechin

Catalog No: #S2856


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Epicatechin
Brief Description	Inhibitors
Purification	94.71%
Target Name	Others
Calculated MW	290.28
Formulation	C15H14O6
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: C1[C@H]([C@H](OC2=CC(=CC(=C2)O)O)C3=CC(=C(C=C3)O)O)O
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Kiewning D, et al. Food Chem. 2013 Sep 1;140(1-2):99-104.

Note: This product is for in vitro research use only