

Vanillyl Alcohol

Catalog No: #S2860

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vanillyl Alcohol
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	154.16
Formulation	C8H10O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: COC1=C(O)C=CC(CO)=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Lahtinen M, et al. *Org Biomol Chem*. 2013 Sep 7;11(33):5454-64.

Note: This product is for in vitro research use only