

## Tyrosol


Catalog No: #S2861

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tyrosol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	138.07
Formulation	C8H10O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: OCCC1=CC=C(O)C=C1  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Chandramohan R, Pari L, Rathinam A, Sheikh BA. Chem Biol Interact. 2015 Mar 5;229:44-54.

Note: This product is for in vitro research use only