

p-Coumaric acid


Catalog No: #S2863

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | p-Coumaric acid |
| Brief Description | Inhibitors |
| Purification | 99.00% |
| Target Name | Others |
| Calculated MW | 164.16 |
| Formulation | C ₉ H ₈ O ₃ |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Other SMILES: C1=CC(=CC=C1/C=C/C(=O)O)O Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Dang YJ, Zhu CY. J Food Sci. 2015 Mar;80(3):C504-9.

Note: This product is for in vitro research use only