

## Oleanolic Acid

Catalog No: #S2865


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Oleanolic Acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	456.72
Formulation	C30H48O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OtherSMILES:C[C@]12CC[C@@H](C[C@@H]1CC[C@@]3[C@@H]2CC=C4[C@]3(CC[C@@]5[C@H]4CC(CC5)(C)C)C(=O)O)C)(C)CO

Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Ethanol: 7 mg/mL (15.32 mM); DMSO: 21 mg/mL (45.98 mM)Appearance:White

PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Liu J, J Ethnopharmacol, 1995, 49(2), 57-68.

Note: This product is for in vitro research use only