

Echinocystic acid

Catalog No: #S2866


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Echinocystic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	472.7
Formulation	C30H48O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: C[C@]12CC[C@@H](C[C@@H]1CC[C@@]3[C@@H]2CC=C4[C@]3(C[C@H]([C@@]5[C@H]4CC(CC5)(C)C(=O)O)O)C)C)O

Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Lai P, Liu Y. Life Sci. 2014 Oct 2;114(2):62-9.

Note: This product is for in vitro research use only