

Diosgenin

Catalog No: #S2867


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Diosgenin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	414.63
Formulation	C27H42O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:C[C@@H]1CC[C@@]2([C@H]([C@H]3[C@@H](O2)C[C@@H]4[C@@]3(CC[C@H]5[C@H]4CC=C6[C@@]5(CC[C@@H](C6)O)C)C)OC1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Ethanol: 83 mg/mL (200.18 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Wu L, et al. Biol Reprod. 2015 Jan;92(1):24.

Note: This product is for in vitro research use only