

Quercitrin

Catalog No: #S2871


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Quercitrin
Brief Description	Inhibitors
Purification	99%+
Target Name	Others
Calculated MW	448.39
Formulation	C21H20O11
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:C[C@H]1[C@@H]([C@H]([C@H]([C@@H](O1)OC2=C(OC3=CC(=CC(=C3C2=O)O)O)C4=CC(=C(C=C4)O)O)O)OPathways:
OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Ma JQ, et al. Food Funct. 2015 Oct 28.

Note: This product is for in vitro research use only