

## Ginsenoside Rg2

Catalog No: #S2873


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Ginsenoside Rg2
Brief Description	Inhibitors
Purification	98%+
Target Name	GSK-3 antagonist
Calculated MW	785.01
Formulation	C42H72O13
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Cancer SMILES: CC1OC(OC2C(O)C(O)C(CO)OC2OC2CC3(C)C(CC(O)C4C(CCC34)C(C)(O)CCC=C(C)C)C3(C)CCC(O)C(C)(C)C23)C(O)C(O)C1O

Pathways: PI3K/Akt/mTOR signaling; Stem Cells Receptor: GSK-3β

Boling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Cho YS, et al. Korean J Physiol Pharmacol. 2013 Apr;17(2):133-7.

Note: This product is for in vitro research use only