

Isopsoralen

Catalog No: #S2875


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Isopsoralen
Brief Description	Inhibitors
Purification	99.98%
Target Name	Others
Calculated MW	186.17
Formulation	C11H6O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:O=C1OC2=C(C=CC3=C2C=CO3)C=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lu H, et al. Afr Health Sci. 2014 Sep;14(3):641-7.

Note: This product is for in vitro research use only