

Flavone


Catalog No: #S2876

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Flavone
Brief Description	Inhibitors
Purification	99.00%
Target Name	NOS antagonist
Calculated MW	222.25
Formulation	C15H10O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: O=C1C=C(OC2=CC=CC=C2)C1=CC=CC=C1
Pathways: Immunology/Inflammation
Receptor: NOS
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Zhu W, et al. Biochem Biophys Res Commun. 2015 Mar 13;458(3):590-5.

Note: This product is for in vitro research use only