

Mucic Acid

Catalog No: #S2877


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mucic Acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	210.14
Formulation	C6H10O8
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:[C@@H]([C@@H]([C@H](C(=O)O)O)O)([C@@H](C(=O)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Mehti? T, et al. Carbohydr Res. 2015 Jan 30;402:102-10.

Note: This product is for in vitro research use only