

Ginsenoside Rd

Catalog No: #S2878


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ginsenoside Rd
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	947.15
Formulation	C48H82O18
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:CC(C)=CCCC(C)(OC1OC(CO)C(O)C1O)C1CCC2(C)C1C(O)CC1C3(C)CCC(O)C(C)(C)C3C(CC21C)OC1OC(CO)C(O)C(O)C1OC1OC(CO)C(O)C(O)C1OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Wang B, et al. Neural Regen Res. 2014 Sep 15;9(18):1678-87.

Note: This product is for in vitro research use only