

Ecdysterone

Catalog No: #S2881


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ecdysterone
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	480.65
Formulation	C27H44O7
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:C[C@]12CC[C@H]3C(=CC(=O)[C@H]4[C@@]3[C[C@@H]([C@@H](C4)O)O)C[C@@]1(CC[C@@H]2[C@](C)([C@@H](CCC(C)(C)O)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Ethanol: 96 mg/mL (199.73 mM); DMSO: 96 mg/mL (199.73 mM)Appearance:Faint Gray PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Kizelsztejn P, et al. Am J Physiol Endocrinol Metab, 2009, 296(3), E433-439.

Note: This product is for in vitro research use only