

Genistin


Catalog No: #S2882

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Genistin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	432.37
Formulation	C ₂₁ H ₂₀ O ₁₀
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: C1=CC(=CC=C1C2=C(C(=O)O)O[C@H]4[C@@H]([C@H]([C@@H]([C@H](O4)CO)O)O)O Pathways: Other
 sReceptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Xiong YJ, et al. Arch Pharm Res. 2013 Mar;36(3):345-52

Note: This product is for in vitro research use only