

Syringic acid


Catalog No: #S2883

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Syringic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	198.18
Formulation	C9H10O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: COC1=CC(=CC(OC)=C1)C(=O)O
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Gven M, et al. Turk J Med Sci. 2015;45(1):233-40.

Note: This product is for in vitro research use only