

Esculin

Catalog No: #S2885


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Esculin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	340.28
Formulation	C15H16O9
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:C1=CC(=O)OC2=CC(=C(C=C21)O[C@H]3[C@@H]([C@H]([C@@H]([C@H](O3)CO)O)O)O)OPathways:OthersReceptor:Others
 Boling pt:Melting pt:Solubility:Ethanol: 13 mg/mL (38.2 mM); Water: 2 mg/mL (5.87 mM); DMSO: 68 mg/mL (199.83 mM)Appearance:White
 PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Mokdad-Bzeouich I, et al. J Antibiot (Tokyo). 2015 Mar;68(3):148-52.

Note: This product is for in vitro research use only