

Aleuritic Acid

Catalog No: #S2886


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Aleuritic Acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	304.43
Formulation	C16H32O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:OCCCCC(O)C(O)CCCCC(O)=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Nyasse B, et al. Nat Prod Res. 2006 Apr;20(4):391-7.

Note: This product is for in vitro research use only