

Trigonelline

Catalog No: #S2887


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Trigonelline
Brief Description	Inhibitors
Purification	99%+
Target Name	Others
Calculated MW	137.14
Formulation	C7H7NO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: C[N+]1=CC=CC(=C1)C(=O)O
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Zhang DF, et al. Asian Pac J Trop Med. 2015 Aug;8(8):651-4.

Note: This product is for in vitro research use only