

Hordenine

Catalog No: #S2890


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Hordenine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	165.23
Formulation	C10H15NO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other SMILES: CN(C)CCC1=CC=C(O)C=C1 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Ethanol: 33 mg/mL (199.72 mM); DMSO: 33 mg/mL (199.72 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Kim SC, et al. Food Chem. 2013 Nov 1;141(1):174-81.

Note: This product is for in vitro research use only