

Monoammonium glycyrrhizinate

Catalog No: #S2891


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Monoammonium glycyrrhizinate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	839.97
Formulation	C42H62O16H3N
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: C[C@]12CC[C@](C[C@H]1C3=CC(=O)[C@@H]4[C@]5(CC[C@@H](C([C@@H]5CC[C@]4([C@@]3(CC2)C)C)(C)C)O[C@@H]6[C@@H]([C@H]([C@@H]([C@H](O6)C(=O)O)O)O)[C@H]7[C@@H]([C@H]([C@@H]([C@H](O7)C(=O)O)O)O)C)(C)C(=O)O.N Pathways: Other
 sReceptor: Others Boling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Shi JR, et al. Int Immunopharmacol. 2010 Oct;10(10):1235-41.

Note: This product is for in vitro research use only