

## Norcantharidin

Catalog No: #S2894


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Norcantharidin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	168.15
Formulation	C <sub>8</sub> H <sub>8</sub> O <sub>4</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: O=C1OC(=O)C2C3CCC(O3)C12  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Yuan X, et al. Pharmacology. 2015;95(1-2):1-9.

Note: This product is for in vitro research use only