

## Helenine

Catalog No: #S2896


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Helenine
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	232.33
Formulation	C15H20O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
 SMILES: CC1=C[C@@H]3[C@@H](C2)OC(=O)C3=C)C  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
 Appearance: White Crystal  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Shope RE. J Exp Med. 1966 Feb 1;123(2):213-27.

Note: This product is for in vitro research use only