

Shikonin

Catalog No: #S2897


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Shikonin
Brief Description	Inhibitors
Purification	98%+
Target Name	Others
Calculated MW	288.29
Formulation	C16H16O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CC(C)=CCC(O)C1=CC(=O)C2=C(O)C=CC(O)=C2C1=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Red PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Zhao Q, et al. Int J Y. 2015 Sep 15;137(6):1446-56.

Note: This product is for in vitro research use only