

## Liquiritin

Catalog No: #S2899


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Liquiritin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	418.39
Formulation	C <sub>21</sub> H <sub>22</sub> O <sub>9</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Other SMILES: C1[C@H](OC2=C(C(=O)C=CC(=C2)O)C3=CC=C(C=C3)O[C@H]4[C@@H]([C@H]([C@@H]([C@H](O4)CO)O)O)O Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Zhang W, et al. Biomed Chromatogr. 2014 Sep;28(9):1271-7.

Note: This product is for in vitro research use only