

Daidzin

Catalog No: #S2901


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Daidzin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Dehydrogenase inhibitor
Calculated MW	416.37
Formulation	C21H20O9
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: C1=CC(=CC=C1C2=C(C(=O)C=CC(=C3O[C@H]4[C@@H]([C@H]([C@@H]([C@H](O4)CO)O)O)O)OP Pathways: Metabolism
 Receptor: ALDH-; ALDH-Boling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Zaheer M, Reddy VD, Giri CC. Nat Prod Res. 2015 Jul 8:1-6.

Note: This product is for in vitro research use only