

Sclareolide

Catalog No: #S2904


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sclareolide
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	250.38
Formulation	C16H26O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:C[C@]12CCCC([C@@H]1CC[C@]3([C@@H]2CC(=O)O3)C)(C)CPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Ethanol: 50 mg/mL warmed (199.69 mM); DMSO: 50 mg/mL (199.69 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Jasinski M, et al. Plant Cell, 2001, 13(5), 1095-1107.

Note: This product is for in vitro research use only