

5alpha-Cholestan-3-one

Catalog No: #S2905


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	5alpha-Cholestan-3-one
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	386.67
Formulation	C27H46O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:C[C@H](CCCC(C)C)[C@H]1CC[C@@H]2[C@@]1(CC[C@H]3[C@H]2CC[C@@H]4[C@@]3(CCC(=O)C4)C)CP

Pathways:Others
 Receptor:Others
 Boling pt:Melting pt:Solubility:DMSO: 10 mM
 Appearance:White Powder
 Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Kasimov MR, et al. Biochim Biophys Acta. 2015 May;1851(5):674-85.

Note: This product is for in vitro research use only