

## Tanshinone IIA

Catalog No: #S2906


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Tanshinone IIA
Brief Description	Inhibitors
Purification	99.50%
Target Name	VEGFR antagonist
Calculated MW	294.35
Formulation	C19H18O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: CC1=COC2=C1C(=O)C(=O)C1=C2C=CC2=C1CCCC2(C)C  
Pathways: Angiogenesis; Tyrosine Kinase/Adaptors Receptor: VEGFR2  
Boiling pt: Melting pt: Solubility: DMSO: 5 mg/mL (16.98 mM)  
Appearance: Red Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Fu J, et al. Eur J Pharmacol, 2007, 568(1-3), 213-221.

Note: This product is for in vitro research use only