

Palmitic acid


Catalog No: #S2908

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Palmitic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	256.42
Formulation	C16H32O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CCCCCCCCCCCCCCCC(O)=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Nagata Y, et al. Leuk Res. 2015 Jun;39(6):638-45.

Note: This product is for in vitro research use only