

Fraxetin

Catalog No: #S2909


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fraxetin
Brief Description	Inhibitors
Purification	99%+
Target Name	Others
Calculated MW	208.17
Formulation	C10H8O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: COC1=C(O)C(O)=C2OC(=O)C=CC2=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Wang H, et al. Mol Med Rep. 2014 Nov;10(5):2341-5.

Note: This product is for in vitro research use only