

Stevioside

Catalog No: #S2911


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Stevioside
Brief Description	Inhibitors
Purification	98.00%
Target Name	TLR inhibitor
Calculated MW	804.88
Formulation	C38H60O18
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology
 SMILES: C[C@@]12CCC[C@@]1([C@H]1CC[C@]34[C@H]2CC[C@](C3)(C(=C)C4)O[C@H]5[C@@H]([C@H]([C@@H]([C@H](O5)CO)O)O)[C@H]6[C@@H]([C@H]([C@@H]([C@H](O6)CO)O)O)(C)C(=O)O[C@H]7[C@@H]([C@H]([C@@H]([C@H](O7)CO)O)O)OP
 Pathways: Immunology/Inflammation
 Receptor: TLR2
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Wang T, et al. Int Immunopharmacol. 2014 Sep;22(1):192-9.

Note: This product is for in vitro research use only