

Ergosterol


Catalog No: #S2912

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ergosterol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	396.66
Formulation	C28H44O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: C[C@H](/C=C/[C@H](C)C(C)C)[C@H]1CC[C@@H]2[C@@]1(CC[C@H]3C2=CC=C4[C@@]3(CC[C@@H](C4)O)C)C Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 0.01 mg/mL (0.02 mM) Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Haubrich BA, et al. J Lipid Res. 2015 Feb;56(2):331-41.

Note: This product is for in vitro research use only